

# CacoGoblet, a ready-to-use intestinal cell-based model to screen compounds' anti-inflammatory capabilities

Jonatan Cucala, Itziar Garzón, Marta Ollé, Gemma Montoya and Lourdes Gombau  
MedTech Barcelona, Barcelona, Spain (contact: [jucucala@medtechbcn.com](mailto:jucucala@medtechbcn.com))

## INTRODUCTION


Identification of active molecules with gastrointestinal anti-inflammatory properties mainly relies on time-consuming methods that quantify the effect of these compounds on the expression and/or release of cytokines (e.g., IL-8). To speed up this process, a 21-day differentiated co-culture of human colon carcinoma cell lines (Caco-2 and HT-29) grown in 24-insert permeable supports was used to screen compound gastrointestinal (GI) anti-inflammatory properties. The approach is based on the correlation between gut inflammation and the evaluation of two indicators of cell barrier integrity, the Transepithelial Electrical Resistance (TEER) and the Lucifer Yellow Paracellular Permeability (LY  $P_{app}$ /Flux).

## OBJECTIVES


The main objective of this study is to evaluate the effectiveness of measuring two indicators of cell barrier integrity (TEER and LY  $P_{app}$ /Flux) in a co-culture of Caco-2 and HT-29 cells to assess compound anti-inflammatory properties.

## MATERIALS AND METHODS

### LIQUEFACTION AND EXCHANGE OF SHIPPING MEDIUM (SM)



### TESTING ANTI-INFLAMMATORY ACTIVITY



## RESULTS

### EVALUATION OF TEER AND LY $P_{app}$ /Flux

Quantifying TEER and LY  $P_{app}$ /Flux in cytokine-induced cells not exposed to baicalein reduced TEER values by 35% and increased LY  $P_{app}$ /Flux by 50%. Conversely, baicalein-treated cells normalized cell barrier integrity in a dose-dependent manner.




Figure 1. TEER measurements in 72h cytokine-induced Caco-2/HT-29 cells incubated in the absence/presence of increasing concentrations of baicalein, a reference anti-inflammatory compound.

### INTERLEUKIN-8 (IL-8) SECRETION

The disruption and recovery of cell barrier integrity was correlated with the release of IL-8. Results indicated that IL-8 release decreased with increasing concentrations of baicalein, correlated with the evolution of TEER and LY  $P_{app}$ /Flux values in a healthy barrier.




Figure 3. Measurement of IL-8 secretion by enzyme-immunoabsorbent assay (ELISA) after 72h of Caco-2/HT-29 exposure to a cocktail of cytokines in the absence/presence of increasing concentrations of baicalein.




Figure 2. LY  $P_{app}$ /Flux in 72h cytokine-induced Caco-2/HT-29 cells incubated in the absence/presence of increasing concentrations of baicalein, a reference anti-inflammatory compound. Bars indicate LY  $P_{app}$  and dots indicate LY Flux values.

### REPRODUCIBILITY BETWEEN ASSAYS




Figure 4. TEER and LY  $P_{app}$  values up to 72h incubation of cytokine-induced cells in the absence/presence of 100μM baicalein. Bars indicate TEER and dots indicate LY  $P_{app}$  values.

Measurement of two indicators of cell barrier integrity (TEER and LY  $P_{app}$ /Flux) in a co-culture of Caco-2 and HT-29 cells can be used to replace IL-8 secretion for screening compounds with anti-inflammatory activity at the early stages of drug development.



#### REFERENCES:

- Putt et al., Food & Function, 8(1), 821-834, 2017.
- Ponce de León-Rodríguez et al., Critical Reviews in Food Science and Nutrition, 59(22), 3648-3666, 2019.
- Vancamelbeke & Vermeire, Expert Review of Gastroenterology & Hepatology, 11(9), 821-834, 2017.
- Joshi et al., In Vitro Models, 1(3), 213-227, 2022.